


System AmbiTherm N4 – Pokyny k zabalení pro GO! Ambiente (15 až 25 °C)



Č.	Popis	Počet	Teplota
1	Krabice 0201w 448x393x318mm	1	22 °C ± 5 °C
2	Meditherm 4.0 Base 15g/litr	1	22 °C ± 5 °C
3	Medicool MC0215 +20 °C	7	20 °C ± 2 °C
4	Krabice 0201B 235x180x95mm	1	22 °C ± 5 °C
5	Zatížení 235x180x95mm	1	20 °C ± 2 °C
6	Meditherm 4.0 Lid 15g/litr	1	22 °C ± 5 °C